


UV Transilluminator



Key Features

- High Quality UV light source
- Electronic Starter UV light tubes to ensure long lasting working capabilities
- Cooling fan cum exhaust
- Orange acrylic UV shield for safety
- Works best with Ethidium Bromide

Specifications

Cat No.	BB-UVT-B	BB-UVT-A	BB-UVT-I	BB-UVT-D
UV Source	4 x 4 Watts	4 x 8 Watts	6 x 8 Watts	2 x 4 x 8 watts
Wavelength	365nm	312/254nm	312/254/365	312/254/365 + White Light
Filter	Indigenous	Indigenous	Imported	Imported + White Light
UV Shield Cover	UV Protection	UV Protection	UV Protection	UV Protection
Dimensions	20 x 10 cm	20 x 20 cm	20 x 20 cm	2x (20x20cm)