


Vertical Gel Electrophoresis Systems(PAGE)



Key Features

- Basic Unit made of high-quality acrylic material
- All units include spirit level, platinum electrodes, glass plates – notched and rectangular (2 set), fixed rubber gasket, safety lid, clamps and screws (2 sets in dual Model), connecting cords
- Combs and spacers made of high quality Teflon
- Available in single and Dual Models
- Water Circulation (only in dual models)

Specifications

Cat No.	BB-ES-V-Mini	BB-ES-V-Midi	BB-ES-V-Maxi
Dimension (LXBXH)	15 x 18 x 15 cm	25 x 20 x 27 cm	29 x 20 x 21 cm
No. Of Samples	7	13	
Gel Size	8 x 7 cm	16 x 14 cm	16 x 20 cm
Buffer Volume	100 + 150 ml	150 + 200 ml	300 + 600 ml
Combs	7 well 0.5mm	13 well 0.5, 1, 1.5mm	20 well 0.5, 1, 1.5mm (2)
Spacers	0.5, 1, 1.5 mm (2)	0.5, 1, 1.5 mm (4)	0.5, 1, 1.5 mm (6)

Cat No.	BB-ES-DV-Mini	BB-ES-DV-Midi	BB-ES-DV-Maxi
Dimension (LXBXH)	16 x 20 x 15cm	23 x 18 x 27 cm	29 x 22 x 23 cm
No. Of Samples	14	26	40
Gel Size	8 x 7 cm (2 Gel)	16 x 14 cm (2 Gel)	16 x 20 cm (2 Gel)
Buffer Volume	150 + 200 ml	250 + 500 ml	300 + 700 ml
Combs	7 well 0.5mm, 1 mm (2)	13 well 0.5, 1mm (2)	20 well 0.5, 1mm (2)
Spacers	0.5, 1 mm(4)	0.5, 1 mm(4)	0.5, 1 mm(4)