

VORTEX MIXER



BB-VS-C



BB-VS-A

Stirring & Mixing

Key Features

- The unit is meant for quick and uniform mixing of liquids in test tubes, small flasks or bottles
- The mixer consists of an enamel finished housing containing the motor with an eccentric drive to a rubber buffer on the top
- Specially designed vacuum suction feet for body stability

Specification

Cat No.	BB-VS-A	BB-VS-C
Power	60 W	60W
Mixing motion	Orbital	
Orbital Diameter	4 mm	4 mm
Motor Type	Shaded Pole	Shaded Pole
Speed Range	Fixed speed: 2500RPM	Adjustable: 0-2800 RPM
Run Type	Touch/ Continuous	Touch/ Continuous
Weight	3.5 Kg	3.5Kg
Accessories	Top Cup (Default), Plate Adapter, Sponge Adapter for Tubes	Top Cup (Default), Plate Adapter, Sponge Adapter for Tubes