

Bottle-top Dispenser



Key Features

- Excellent chemical resistance, components are made of PTFE, FEP, BSG, PP
- Fully autoclavable at 121°C
- Four ranges of bottle-top dispenser covering a volume range from 0.5mL to 50mL
- Easy for cleaning and maintenance
- The optional flexible discharge tube with safety handle permits fast and precise dispensing
- Vapor pressure max. 500mbar, viscosity max. 500mm²/s, temperature max. 40°C, density max. 2.2g/cm³
- Supplied with S40, GL32, GL38, GL25, GL28
- Adapters - (GL-32 – 45/32mm, GL-38 – 45/38mm, GL25- 32/25 mm, GL28 – 32/25 mm, s40-45/40mm)

Specification

Cat No.	Volume Range (ml)	Graduation (ml)	Accuracy (%)	Precision (%)
BB-BTD-V1	0.5-50	0.1	±0.5	0.2
BB-BTD-V2	1.0-10	0.2	±0.5	0.2
BB-BTD-V3	2.5-25	0.5	±0.5	0.2
BB-BTD-V4	5.0-50	1.0	±0.5	0.2