

Electrophoresis Power Supply- Analog



Key Features

- Microprocessor Controlled,
- 6 step Voltage Selectable,
- Running Status Led,
- Buzzer alert on set timer end,
- 2 in Parallel – Output Terminal,
- Android Mobile App (Bluetooth) Controlled ,20 to 25 mts Distance Access
- Future upgrade to 9 Step on Extra Cost. On-site Calibration
- Android Mobile App – On extra Cost

Specifications

Cat No.	BB-PSB1	BB-PSB2	BB-PSB3
Steps	6 Step Voltage Selection		
Voltage	50V, 100V, 150V, 200V, 250V,300V		
Current	100mA	200mA	300mA
Operation	LED Display		
Control Type	Soft key		
Safety	Short Circuit Protection and Overload Protection		
Input voltage	AC 220 V, 50 Hz		