

# WATER BATH



BB-WB-4H



BB-WB-8H

Thermal Controller

## Key Features

- New design structure, which could avoid the leakage of steam
- High precision on controlling the temperature
- digital display with key operating button

## Specifications

Cat No.	BB-WB-4H	BB-WB-6H	BB-WB-8H
<b>Column</b>	2 rows 4 holes	2 rows 6 holes	2 rows 8 Holes
<b>Power rating</b>	800W	1200W	1400W
<b>Inner Chamber size (W x L x H) mm</b>	305 X 305 X 130	305 X 470 X 130	305 x 620 x 130
<b>Internal volume (approx.)</b>	12 L	18 L	24L
<b>Temp. range</b>	Rt+5~100°C		
<b>Temp. Resolution Ratio</b>	0.1°C		
<b>Temp. Accuracy</b>	±0.5°C		
<b>Water sink</b>	Stainless steel		
<b>Outer shell</b>	surface electrolysis treatment Cold rolled steel plate		
<b>Sensor</b>	Thermal Resistor		
<b>Safety device</b>	Over-temperature alarm, short-circuit protection		