


Ph Meter



Key Features

- Microprocessor Controlled
- Jumbo – LCD display
- ABS – Plastic Body
- Android Mobile App Control
- Bluetooth Distance – 20 to 25mts
- Soft Key & Buzzer Key Alert
- Inbuilt Magnetic Stirrer Control
- No external power require for Stirrer
- Stirrer Speed 0 to 1200 rpm / 510 step
- Below 0 pH & Alert Alarm LED
- On-site Calibration
- Calibration is retained in the memory even after the instrument is switched off

Specifications

Cat No.	BB-PH-3P	BB-PH-5P	BB-PH-8P
Measure Range	0 to 14 pH		
Resolution	0.01 pH		
Resolution	0.1 mV		
Resolution	.1°C		
Accuracy and Repeatability	± 0.03		
mV Range	+ 414.1 to - 414.1		
Accuracy	± 0.1 C		
Calibration	3 Point	5 Point	8 Point
Measures	pH & mV, Temperature		
Temperature Compensation	Auto/Manual, ATC (PT100 Sensor – 3 Wire)		
Auto / Manual Temp Range	0 °C to 99.9°C		
Input voltage	AC 220 V / 50 Hz		