

# Rotary Shaker



## Key Features

- Speed: 40 to 250 rpm
- Control Type: PID – Rpm Control
- Top platform: Square Type
- Cabinet: GI body with Powder Coated
- Top View: GI top view
- Platform: Flask Holder
- Slow start, slow stop
- Additional flask holder 500ml, 1000ml
- Input voltage: AC 220 V, 50 Hz

## Specification

Cat No.	BB-RS-A	BB-RS-D
RPM	40-250 RPM	
Accuracy	±1 RPM	
Flask Holder	Universal/Clamps (16x250 ml)	
Flask Holding platform	380 x 380 mm	
Timer	1 min – 99hours 59 min	
Control	Knob Control	PID control with Digital Display