

HORIZONTAL AUTOCLAVE

Key Features

- Triple walled with separate steam jacket and boiler
- Inner chamber and steam jacket are made of heavy gauge Stainless Steel sheet with leak proof argon arc welding
- The sterilizer has single piece door made of thick stainless Steel plate.
- All the sterilizers are hydraulically tested to withstand pressure upto 40 psi.
- Mounted on tubular steel frame with ground leveling screwed flanges.
- The space between outer wall and steam jacket is filled with glass wool to minimize the heat losses
- Made of heavy gauge Stainless Steel Sheet
- Fitted with flanged immersion type heaters at the back
- Also fitted with Low Water Level Cut off Device to protect the heaters from burning out dry
- Fitted with Water Level Indicator, Water Inlet and Outlet Valves
- Includes Automatic pressure cut off switch, digital temperature indicator cum controller, Vacuum Breaker, Water level indicator, moisture trap and multiport valve



Specifications

Cat No.	Inner Dimensions (Dia X Depth)	Capacity (Approx.)
BB-HA-78	400 x 600mm – 16" X 24"	78 Ltrs
BB-HA-115	400 x 900mm – 16" X 36"	115 Ltrs
BB-HA-162	400 x 1200mm – 16" X 48"	162 Ltrs
BB-HA-221	500 x 1200mm – 20" X 48"	221 Ltrs
BB-HA-300	550 x 1100mm – 22" X 48"	300 Ltrs
BB-HA-350	600 x 1200mm – 24" X 48"	350 Ltrs