

MicroPipette A1



Key Features

- Fully Autoclavable
- Comprehensive Volume range from 0.1µL to 10 ml
- Available as both fixed and variable volume
- 8 and 12 channel Pipette options
- Seamlessly rotating dispensing heads for effortless pipetting
- Spring loaded tip cones for easy cleaning and maintenance

Specifications

VARIABLE VOLUME A1 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A1-V1	0.1-2.5	0.05	±3	3.00
BB-A1-V2	0.5-10	0.1	±1.5	1.50
BB-A1-V3	2-20	0.5	±1.2	1.00
BB-A1-V4	5-50	0.5	±0.9	0.60
BB-A1-V5	10-100	1	±1	0.
BB-A1-V6	20-200	1	±0.8	0.3
BB-A1-V7	100-1000	5	±0	
BB-A1-V8	1000-5000	50	±0.6	0.30

FIXED VOLUME A1 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A1-F1	5	-	±1.3	1.20
BB-A1-F2	10	-	±0.8	0.80
BB-A1-F3	20	-	±0.6	0.50
BB-A1-F4	50	-	±0.5	0.3
BB-A1-F5	100	-	±0.5	0.30
BB-A1-F6	200	-	±0.4	0.20
BB-A1-F7	500	-	±0.3	0.2
BB-A1-F8	1000	-	±0.3	0.20
BB-A1-F9	2000	-	±0.3	0.15
BB-A1-F10	5000	-	±0.3	0.15

8 CHANNEL A1 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A1-8C1	0.5-10	0.1	±2.50	2.5
BB-A1-8C2	5-50	0.5	±1.50	1.0
BB-A1-8C3	50-300	5	±1.00	0.

12 CHANNEL A1 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A1-12C1	0.5-10	0.1	±2.50	2.50
BB-A1-12C2	5-50	0.5	±1.50	1.00
BB-A1-12C3	50-300	5	±1.00	0.5

MicroPipette A2



Key Features

- Semi Autoclavable (Tip Cone & Tip Ejector)
- Comprehensive Volume range from 0.1µL to 10 ml
- Available as both fixed and variable volume
- 8 and 12 channel Pipette options
- Seamlessly rotating dispensing heads for effortless pipetting
- Spring loaded tip cones for easy cleaning and maintenance

Specifications

VARIABLE VOLUME A2 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A2-V1	0.1-2.5	0.05	±3	3.00
BB-A2-V2	0.5-10	0.1	±1.5	1.50
BB-A2-V3	2-20	0.5	±1.2	1.00
BB-A2-V4	5-50	0.5	±0.9	0.60
BB-A2-V5	10-100	1	±1	0.40
BB-A2-V6	20-200	1	±0.8	0.
BB-A2-V7	100-100			
BB-A2-V8	1000-5000	50	±0.6	0.30

FIXED VOLUME A2 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A2-F1	5	-	±1.3	1.20
BB-A2-F2	10	-	±0.8	0.80
BB-A2-F3	20	-	±0.6	0.50
BB-A2-F4	50	-	±0.5	0.3
BB-A2-F5	100	-	±0.5	0.30
BB-A2-F6	200	-	±0.4	0.2
BB-A2-F7	500	-	±0.3	0.20
BB-A2-F8	1000	-	±0.3	0.
BB-A2-F9	2000	-	±0.3	0.15
BB-A2-F10	5000	-	±0.3	0.1

8 CHANNEL A2 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A2-8C1	0.5-10	0.1	±2.50	2.5
BB-A2-8C2	5-50	0.5	±1.50	1.0
BB-A2-8C3	50-300	5	±1.00	0.5

12 CHANNEL A2 PIPETTE				
Cat No.	Volume Range (µl)	Increment (µl)	Accuracy (%)	Precision (%)
BB-A2-12C1	0.5-10	0.1	±2.50	2.50
BB-A2-12C2	5-50	0.5	±1.50	1.00
BB-A2-12C3	50-300	5	±1.00	0