

Pipette Controller



Key Features

- Compatible with most of the serological and glass pipettes from 0.1 -100mL
- Six speeds (only for + Model) of selection for aspiration and dispensing different liquids
- Large LCD display (only for + Model) showing low battery warning and speed settings
- Enables single-handed operation with minimum effort
- Light and ergonomic design providing easy usability
- High capacity Li-ion battery enable long runtime of operation
- Quiet & Powerful pump fills a 25mL pipette in <5 seconds
- 0.45µm replaceable hydrophobic filter
- Rechargeable during use

Specification

Cat No.	BB-PG+	BB-PG
Aspiration Speed	25mL<7s	25mL<6s
Dispensing Speed	Motor (6 shift)/Gravity	Motor/Gravity
Battery	Lithium-ion	
Battery Service Life	More than 8 Hours Intermittence Use	
Charging Time	2-3 Hours	
Pipette Types	Serological or Plastic Pipette (0.1-100mL)	
Filter	0.45µm Hydrophobic	
Weight	200g	