

# Fume Hood



## Key Features

- Fabricated out of Heavy gauge MS duly powder coated
- These are designed to throw-out all toxic/harmful vapours thus protecting costly instruments and the person working in the lab
- Working table top is made of stainless steel and a small wash basin with connections provided for inlet and outlet of water
- The front door moves vertically up and down with concealed counter balanced weight
- The unit is fitted with fluorescent light and a gas cock for gas/air supply
- The front facing panel is fitted with 15/5 Amp socket with switches for exhaust system and fluorescent light
- To work on 220 / 230 volts AC Supply

## Specifications

Cat No.	BB-FM-3	BB-FM-4	BB-FM-6
<b>Workspace (W x D x H)</b>	(2'x2'x2')	(3'x2'x2')	(4'x2'x2')
<b>Airflow Pattern</b>	100% Exhaust		
<b>Noise</b>	<65dB		
<b>Power Supply</b>	220V/50Hz		
<b>Optional Accessories</b>	Front Cover mounting base with storage space and door		