

Vacuum Pump



Key Features

- Our vacuum pump is a piston-powered, oil-free pump. With innovative electronic, mechanical technology and human design concept, compact and light weight, clean and maintenance free, safe and comfortable
- Most ideal for use with Rotary Vacuum Evaporators, Vacuum Ovens for continuous operations without the fear of sample contaminations through oil vapours
- Useful for general laboratory applications where even a slightest environment/sample contamination through oil mist is forbidden
- 100% Oil Less Vacuum Pump

Specifications

Cat No.	BB-VP-17	BB-VP-34
Flow rate [L/min] at atm. pressure	17 L/Min	34 L/Min
Ultimate vacuum [mbar abs.]	100 mbar	
Max. Power P [W]	35W	60W
Max. current [A]	0.3	0.4
Motor speed [rpm]	1450	
Hose connections [mm]	9	
Noise [dB]	50	