

Electrophoresis Power Supply Basic-Digital



Key Features

- Microprocessor with Knob Controlled
- LCD Display
- Running / Mode – Status Led
- Buzzer alert on Error / Timer
- 2 in Parallel Output Terminal
- Auto restart after an AC power failure
- Digital Timer : 1 – 99 Min
- Process complete buzzer alert

Specifications

Cat No.	BB-PSDB1	BB-PSDB2	BB-PSDB3
Steps	9 Step Voltage Selection		
Voltage	50V, 75V, 100V, 125V, 150V, 175V, 200V, 225V, 250V		
Current	100mA	200mA	300mA
Operation	LCD Display		
Control Type	Soft key		
Timer	1 -99 Min		
Cabinet	ABS plastic Body		
Security	Short Circuit, Overload Protection		
Input voltage	AC 220 V, 50 Hz		