

PCR Workstation



Clean Air/Sterilization

Key Features

- Antimicrobial coated stainless steel reduces bacterial growth
- A safety shut-off switch automatically turns the ultraviolet light off when the door is open
- Built-in power outlets allow operation of equipment for use in PCR applications
- Workstation includes a 30 minute timer to control UV decontamination
- Built in Pipette holder

Specifications

Cat No.	BB-PCRW
UV Source	254 nm x 2
Dimensions (BxWxH)	24" x 24" x 30"
Timer	0-60 Minutes-
Control	Digital display with timer control