


Submarine (Horizontal) Gel Electrophoresis



Key Features

- Completely moulded system made of transparent high quality Acrylic material
- All the units include basic unit, platinum electrodes, connecting cords, safety lid, gel running tray, gel casting tray and different sizes of combs

Specifications

Cat No.	BB-ES-H-Mini	BB-ES-H-Midi	BB-ES-H-Maxi
Gel Size	10 x 7 cm	13 x 13 cm	20 x 20 cm
Gel Trays	10x7 cm (1), 5x7 cm (2)	13x13 cm (1), 13x6cm (2), 6.5x6cm (4)	20 x 20 cm (1)
Combs	8 well x 1.5mm (2) 3 well x 3mm gel prep	13 Well X 1.5 mm (1) 8 well x 1.5mm (4) 3 well x 3mm gel prep	26 well x 1.5mm (1) 3 well x 3mm gel prep 3 well x 6mm gel prep
No. Of Samples	8	13	
Gel Casting Tray	Universal		
Dimensions (LXBXH)	20x10x8 cm	27 X 16 X 8 cm	35 x 22 x 8 cm