

TISSUE CULTURE RACK



Key Features

- Made of tubular pipe duly epoxy powder coated
- Each rack of size 48" x 18" is covered with thick unbreakable Plexi Glass / Hylem Sheet
- Distance between each rack is 15"
- Each rack is fitted with florescent tubes with Electronic Boosters for minimum flickering effect
- Control panel with timer is provided in the advanced models
- Rack is fitted on Castor Wheels for easy mobility

Specifications

Cat No.	BB-PT-4	BB-PT-4C	BB-PT-5	BB-PT-5C
Racks	4		5	
Fluorescent tubes	12		16	
UV Tubes	0	3	0	4
Timer	No	Yes	No	Yes