


Electroblotter System (Transfer Unit)

Wet Blotters

Key Features

- Basic Unit made of high-quality acrylic material
- All units include Gel cassettes, Filter Pads, Cushion Pads, Platinum electrodes, connecting cords, safety Lids
- Heat Exchanger (only for dual Models)
- Available in single and Dual Models



Specifications

| Cat No. | BB-ET- Mini | BB-ET-D-Mini | BB-ET-Midi | BB-ET-D-Midi |
|-------------------------|-------------|--------------|------------|----------------|
| Gel Size | 8x7 ml | 8x7 ml (2) | 16 x 14 cm | 16 x 14 cm (2) |
| Buffer Reservoir | >350 ml | >600 ml | >1000 ml | >1500 ml |
| Gel Cassette | 1 | 2 | 1 | 2 |
| Filter Pad | 6 | 6 | 6 | 6 |
| Cushion Pad | 12 | 12 | 12 | 12 |
| Heat Exchanger | NO | YES | NO | YES |

Semi dry Blotters

Key Features

- All units include graphite electrodes, non conducting sheets and connecting cords
- Heat Exchanger (only for dual Models)



Specifications

| Cat No. | BB-ETSD- Mini | BB-ETSD-MiDi |
|-----------------------------|---------------|--------------|
| Gel Size | 10x10 cm | 20x20 cm |
| Electrode Assembly | Graphite 1 | Graphite 1 |
| Non-Conducting Sheet | 1 | 2 |