

# Incubator



Thermal Controller

## Key Features

- Microprocessor controlled Ambient +5°C to 60°C (PT 100 Sensor)
- Digital 4-digit 0-9999 minutes Timer
- PID controller – display SV & PV value
- Temperature resolution of 0.1°C, with accuracy of +/- 0.1 °C (heating) making it possible with two air recirculation fan
- PUF insulated for minimal thermal loss
- Auto Restart (On Power Failure)
- Long life & Low Maintenance
- Energy efficiency (SSR Relay Control)
- Safety Fuse & Short Circuit Protection
- Inner body made up of SS (Grade 304) and Outer body made of GI duly powder coated

## Specifications

Cat No.	Inner Chamber Size	No of Shelves
BB-BI-2	055 LTR (2 CFT)	1
BB-BI-4	112 LTR (4 CFT)	2
BB-BI-6	171 LTR (6CFT)	3
BB-BI-8	226 LTR (8CFT)	3
BB-BI-10	275 LTR (10CFT)	4
BB-BI-12	336 LTR (12CFT)	4

## Optional Accessories

- LCD Programmable RTC – Timer (App Control) on additional cost
- Digital Datalogger on additional cost.