


Gel Rocker



BB-GR-A



BB-GR-D

Key Features

- Top platform: Square Type
- Cabinet: GI body with Powder Coated
- Top View: Plate with Anti slip mat
- Input voltage: 12 V / 2A DC – Adapter
- Platform: 240 mm X 240 mm
- Noise & Vibration less DC Motor,
- Auto Restart (On Power Failure)
- Continuous Duty
- Short Circuit & Overload Protection
- Double Decker option also available

Specification

Cat No.	BB-GR-A	BB-GR-D
Speed	0-60 Oscillations	
Top	Plate with Anti-slip Mat	
Platform	380 x 380 mm	
Timer	-	1 min – 99
Control	Knob Control	PID control with Digital Display