


DRY BATH



Thermal Controller

Key Features

- Microprocessor Controlled,
- Inbuilt temperature sensor concept
- Tube size: 1.5 ml / 2.0 ml X 24 Place
- Single Aluminum Block
- +/- 0.1 Deg C Accuracy
- Short Circuit & Heater Protection



Specifications

Cat No.	BB-DB-01
Tube size	1.5 ml / 2.0 ml / 0.2 ml / 0.5ml
Display	Digital Seven Segment
Sensor	PT100 (3 wire)
Controller	PID Temp Controller
Control	Soft Keys
Temp Range	Ambient +/- 5.0 Deg C to 100.0 Deg C
Outer Body	Gi Powder Coated
Inner Body	Aluminum Block
Input voltage	AC 220 V, 50 Hz

Optional Accessories

Customised aluminium block