

Pipette Aid



Pipette Pump



Key Features

- Widest Volume range from 0.1ml to 100 ml
- Reduce hands-on time
- 3 µm replaceable hydrophobic filter
- Easy to clean and maintain
- Silicon adapter for 0.1 – 100 ml pipette tip
- Volume capacity 2ml, 10ml, 25ml
- Thumb wheel guarantees precise operation
- Resistant to acids and bases
- Colour coded by volume

Specification

Cat No.	Volume
BB-PA-100	Pipette Aid 0.1-100 ml
BB-PP-2	Pipette Pump 2 ml
BB-PP-10	Pipette Pump 10 ml
BB-PP-25	Pipette Pump 25 ml

Stand for Micropipette

Cat No.	Description
BB-MPSL	Linear Stand for 6 Pipettes
BB-MPSC	Round stand for 6 Pipettes

